

FMRI – DIFFICULTIES IN MIND READING

- ◎ fMRI: functional (nuclear) magnetic resonance imaging
- ◎ Neuroimaging: get the structure of the brain
 - Want to know how it works: connection brain parts and brain functions
 - Aim: measure the local „thinking activity”
- ◎ Usage (criticism):
 - Lie detector
 - Neural- and psychological-modell checking (think/know experiment)

How we use it

- ⦿ Human attempts frequent
- ⦿ Well-planned tasks or questions
- ⦿ Measures: order of minutes
 - One measure: order of 5 seconds
 - Measure with and without tasks or question, further investigation based on the difference
- ⦿ Overlap the intensity map and brain image

The basics of fMRI



- ⦿ MRI: interaction between spins and magnetic field
- ⦿ QM based phenomena
 - Classical view is almost satisfactory
- ⦿ Find a „think-activity“-sensitive MRI measureable quantity, measure it, and then reflect to think activity

Get the signal - BOLD

- ⦿ Blood-oxygen-level dependence
- ⦿ Hemoglobin: Fe^{2+} can absorb O_2
 - Hemoglobin + O_2 : diamagnetic molecule
 - Hemoglobin – O_2 : paramagnetic ($S=2$)
- ⦿ Measure the oxygen-flow differences in vein
- ⦿ Determine the connection between BOLD signal and thoughts

Get the signal – O_2 flow

- ⦿ The oxygen flow depends on the communicating intensity
- ⦿ Communication needs energy
 - Neuron cells don't have repository
 - Increased activity needs more energy
 - BOLD signal decreasing and then increasing



Time, accuracy and resolution

- ⦿ Time: depends on the BOLD signal, about less than half a minute
 - Do experiments with the same patient, same time
- ⦿ Accuracy: easily detect maximum of BOLD
- ⦿ Space: resolution: in order of *mm×mm×mm*
 - Problem: the motion of the patient

Increasing space resolution

- No new information with increasing resolution
 - Signal comes from multiple capillars
- Solution: BOLD signal minimum: more localised
 - Longer time or
 - Bigger magnetic field
- No (?) news: cerebral tasks are not well localised

Temporal resolution

- ⊙ $\Delta x \cdot \Delta y \cdot \Delta z \cdot \Delta t \approx \text{const}$
- ⊙ Increasing time resolution does count
- ⊙ (more details) **usable image**
- ⊙ Two tricks to improve:
 - Spin echo – gradient echo
 - EPI: echo planar imaging

Measurement: EOM

- First: apply $\mathbf{B}_0 = (0, 0, B_{z,0}) \Rightarrow \mathbf{M} = (0, 0, M_0)$
- Short RF pulse: $\mathbf{B}(t) = B_1 (\sin(\omega t), \cos(\omega t), 0)$

- $$\frac{d\mathbf{M}}{dt} = \gamma (\mathbf{M} \times \mathbf{B}) + (\mathbf{M}_0 - \mathbf{M}) \tau$$

- $$\frac{dM_z}{dt} = \gamma (\mathbf{M} \times \mathbf{B})_z + \frac{M_0 - M_z}{T_1}$$

- $$\frac{dM_i}{dt} = \gamma (\mathbf{M} \times \mathbf{B})_i - M_i / T_2 \quad i \in \{x, y\}$$

Measurement: relaxation and decoherence

- Solutions after


$$M_z = M_0 \left(1 - e^{-t/T_1}\right)$$

„excitation“:


$$M_i = M_0 \cdot e^{i\gamma B_0 t - t/T_2} = M_0 \cdot e^{i\omega_L t - t/T_2}$$

- We can measure M_i
- Important: $T_1 \neq T_2$ $T_2 > T_1$
- New rotating cylindrical coordinates with frequency: $\omega_L = \gamma B_0$


Measurement: spin-echo




Measurement: spin-echo



Measurement: spin-echo



Measurement: slice-excitation




$$\begin{aligned}\frac{dM_i}{dt}(z,t) &= \\ &= \gamma \left(\mathbf{M} \times \mathbf{B}^*(t) \right)_i \cdot e^{i\gamma(B_0 + G_z \cdot z)} \\ &= \gamma M_0 B_1^*(t) e^{i\gamma B_0 - i\gamma G_z \cdot zt}\end{aligned}$$


$$\begin{aligned}\frac{dM_r}{dt}(z,t) &= \\ &= \gamma M_0 \cdot B_1^*(t) \cdot e^{-i\gamma G_z \cdot z \cdot t}\end{aligned}$$

Measurement: slice-excitation and spin-echo

- Find $B_1^*(t)$ for the desired excitation-distribution: $B_1^*(t) = B_1 \sin(x) / x$ $x = \gamma G_z d \cdot t$




Measurement: another way, the gradient echo




Measurement: EPI

- Echo Planar Imaging: increase time resolution even more!




Pulse sequence




Sampling path in k space

What we get: images



An experiment: caffeine

● caffeine as a contrast booster



Limits

- We cannot answer any why question – only answer the question where
- The brain is 3D – the image is 2D, so more experiments needed
- We can measure only the neurons firing, but not the real activity (block or stimulate?)
- High degree of cerebral plasticity: places may vary

Mind reading

- The map of optic nerves to visual cortex is so „localised”, „continuous”
- Training set: known videos and those fMRI signal
- Unknown video: fMRI signal → video images
- Ability to guess what patient think
- Numerical problem
- Further aim: movement of implants
- Success: camera-eye

Presented movie



Reconstructed movies (AHP)



Highest posterior movies



S1



S2



S3